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**World Innovation**  
Disinfection, ultrasound and drying, all in one single appliance

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The sterilization process is a complex group of well-defined protocols, which must be understood and applied carefully and rigorously. The final result of a sterilization process is strictly dependant on the methods of preparation and packaging of the materials, the efficiency of the apparatus and its correct use.

Therefore, it is important to have the necessary resources and tools available, and it is also important, above all, that all the professionals involved in the process are aware of their determining role, based on the correct application of the procedures.

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The sterilization protocol is set out as follows:

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The various phases require the presence of the operator, as well as handling and transfer of the materials and instruments to be treated.

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# The first multifunction appliance

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## THE FIRST MULTIFUNCTION SYSTEM FOR THE MANAGEMENT OF STERILIZATION PROTOCOLS HAS ARRIVED!

Tecno-Gaz S.p.A., the leader in the world of sterilization, has designed a **new and revolutionary automatic appliance**, able to implement the main phases of the sterilization cycle, thus ensuring absolute compliance with procedures. For the first time, all the priority phases are carried out without human intervention, **without any professional risk and greatly reducing the overall time required**. An appliance which is unique in the world and indispensable in any facility.

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**Very simple-to-use, just select the appropriate cycle and activate the sequence.**

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## DISINFECTION

Once the instruments have been placed in the tub and the lid has been closed, the materials remain static.

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## CLEANING

The equipment automatically starts the ultrasound cleaning phase for a pre-set time.

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## RINSE

After the cleaning phase, the liquid empties from the tub and the rinse phase, using pressure water jets, starts.

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## DRYING

After the rinse phase comes the drying phase, using a mechanism which operates by compressed warm air.

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
## END OF THE PROCESS


The equipment emits a sound signal at the end of the cycle. The instruments are ready for packaging.

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





## BEFORE MULTISTERIL

1. insertion in **disinfection** tub
2. Manual **rinse**
3. Manual transfer to the **ultrasound** tub
4. Manual **rinse**
5. Manual transfer to the **drying system** or manual drying
6. Manual **packaging**
7. **Traceability**
8. **Sterilization**

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## WITH MULTISTERIL

1. Insertion in Multisteril
2. Manual packaging
3. Traceability
4. Sterilization

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# MultiSteril: clear and apparent advantages

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## NO HUMAN INTERVENTION

All the disinfection, cleaning, rinsing and drying phases are carried out automatically without any human intervention, allowing the protocols to be complied with correctly, thus avoiding negligence or operational errors.

## ZERO CONTAMINATION RISK

As the infected materials are not handled or transferred manually, accidents in the workplace are avoided and all processes are managed automatically by Multisteril, using preset sequences which do not involve the operator in any way.

## TIME-SAVING

Today, all the sterilization cycle phases are carried out manually by the operator and implementing all the protocols takes up a lot of their time. With Multisteril, there is no need for the presence of an operator and each phase is automatic.

## GREAT COST-SAVING

The whole process excludes the continuous presence of an operator, the management of the liquids used is automatic and, at the end of the cycle, no cleaning or disinfection is required, leading to great savings in time and overall costs. The product pays for itself in a short period of time.

## PERFECT MANAGEMENT OF ALL THE PHASES

You can be absolutely certain that the protocol is implemented correctly and in compliance with the reference procedures.

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
Compact design, ideal for every sterilization room.

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
### MULTISTERIL OPERATES WITH CONNECTIONS TO THE POWER AND WATER SUPPLIES, AND DOES NOT NEED OTHER CONNECTIONS.

An appropriate container is housed within the apparatus for loading ultrasound liquid, which is then managed automatically by the system. **Multisteril is equipped with an instruments basket.** The tub has been designed to hold normal and surgical trays. **The containment tub holds 6 litres and the lid incorporates both the rinse and drying systems.** Sophisticated electronics control and manage the equipment's functions and offers a selection of 3 programmes:

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#### STANDARD PROGRAMME

- Programme with **20 minutes ultrasound treatment**
- Liquid **concentration 2%**
- **Temperature 60°C**

#### INTENSIVE PROGRAMME

- Programme with **20 minutes ultrasound treatment**
- Liquid **concentration 4%**
- **Temperature 30°C**

#### PROGRAMME 3

- **Free selection** programme

Standard

20min

60°C

2%

Intensivo | Intensive

20min

30°C

4%

Programma 3 | Programme 3

LIBERO / FREE

# MultiSteril: le funzioni

# MultiSteril: its functions

## Extra info

### STANDARD

Cestello / Mesh basket  
Art. SVMA069

Cestello portafrese  
Basket for handpieces  
Art. SVMA072

### OPTIONAL

Supporto tray + 2 tray  
Tray holder + 2 trays  
Art. 2210S

## Product info

Dimensioni [LxAxP]   Dimensions [WxHxD]	v415 x 410 x 450 mm
Peso   Weight	10 kg
Assorbimento max.   Max absorption	1,5 kW
Tensione   Main voltage	230 V - 50 Hz
<b>DIMENSIONI VASCA   TUB DIMENSIONS</b>	
Dimensioni [LxAxP]   Dimensions [WxHxD]	332x120x227 mm