

BUILT-IN SET
Ultrasonic scaler
Ablatore a ultrasuoni
Ultraschall-Zahnsteinentferner
Détartreur à ultrasons
Detartador ultrasónico

Operating instructions
Istruzioni per l'uso
Betriebsanleitungen
Mode d'emploi
Modo de empleo

REF 462.9x.03


ACCESSORIES ERSATZTEILE	AND ACCESSOIRES	SPARE PARTS ET	ACCESSORI E RICHANGI	RICAMBI / ACCESORIOS	ZUBEHÖR UND RECAMBIOS
REF 001.18.03	O-rings / Guarnizioni / O-Ringe / Joints toriques / Guarniciones				
REF 480.01	Insert # AX-S / Inserto # AX-S / Einsatz # AX-S / Insert # AX-S / Inserto # AX-S				
REF 480.02	Insert # BX-S / Inserto # BX-S / Einsatz # BX-S / Insert # BX-S / Inserto # BX-S				
REF 480.04	Insert # PX-S / Inserto # PX-S / Einsatz # PX-S / Insert # PX-S / Inserto # PX-S				
REF 109.00	DYNO Dynamometric wrench / Chiavetta dinamometrica / Drehmomentschlüssel / Clavette dynamométrique / Chaveta dinamométrica				
REF 335.1x	Silicone hose / Tubo siliconato / Silikon Schlauch / Tuyau en silicone / Manguera de silicon				

Satelec® is a registered trademark of Electron Medical System SA, Switzerland.

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DESCRIPTION

Ultrasonic piezo-electric scaler for professional dental use.
Intended use: the ultrasonic scaler can be used for scaling, periodontics, prosthetics and endodontics works employing the TKD inserts listed below or the inserts manufactured by SATELEC®.
This medical device meets the requirements of the European Directives 93/42/EEC (Class IIa) and 2007/47/EC.

WARNINGS

- Before using the device, read the following warnings carefully:
- The device must be used according to these operating instructions.
 - The device must be used by competent and qualified personnel only.
 - The device is not suitable for use with in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide.
 - The ultrasonic scaler must be used exclusively with the inserts supplied or optional, with original spare parts or with inserts approved by the manufacturer.
 - The device must be used with biocompatible inserts according to ISO 10993-1 Standard.
 - Do not change the shape of the inserts by bending or filing them.
 - Do not use worn or damaged inserts.
 - LED radiation: do not stare into beam or view directly with optical instruments.

CONTENTS

TITANUS® SLED ultrasonic scaler, insert # AX-S, insert # BX-S and insert # PX-S (REF 462.00.03), DYNO dynamometric wrench, spare O-ring for the hose connector, silicone hose, USC60 electronic board.

ILLUMINATION

The device is equipped with a highly-efficient LED light source.
In order to replace the lighting element, unscrew the front part of the instrument and remove the tapered lighting element. When positioning the new lighting element, be careful to respect the polarity. If the polarity is inverted, the light does not turn on.

CONNECTION OF THE HANDPIECE

The TITANUS® SLED handpiece (REF 462.00.03) must be used together with the specific silicone hose (REF 335.1x) and with the specific USC60 electronic board (REF 387.00).
After positioning the socket, the handpiece can be fitted to the silicone hose by means of a simple snap connection (Fig.1). Do not attempt to screw.

INSTALLATION

To correctly install and connect the USC60 electronic board and the silicone hose to the dental Unit, please refer to the corresponding Application Notes of the SET.
The power supply unit must meet the requirements of the CEI EN 60601-1 standard.
On installation, the handpiece and the inserts must be sterilized before their first use.

COOLING

Operating water pressure for cooling spray can be 0.6 ÷ 2.0 bar. Spray flow rate is higher than 50 ml/min. Cooling spray must be well-filtered.

As inserts heat up when used dry causing heat damage on the teeth, avoid to work dry or, in any case, use inserts without cooling for a few seconds only.

DESCRIPTION OF VIBRATIONS

All ultrasonic vibrations generated by the TITANUS® SLED transducer move longitudinally along the axis of the insert. Inserts therefore vibrate on a single plane and lateral vibrations, usually common to other instruments, do not exist.
For this reason, during procedure, always hold the insert tangentially to the tooth surface.

Do not direct the insert frontally to the enamel.
Never direct the tip of the insert vertically to the tooth.

FITTING OF INSERTS

Select the right shape of insert according to the treatment to be performed and screw it manually on to the handpiece (Fig.2) by means of the dynamometric wrench.
Do not fit inserts while the scaler is on.

MAINTENANCE OF THE INSERTS

The shape and weight of each insert are the result of an exact design and manufacturing so as to obtain maximum performance from the ultra-sonic transducer. Therefore any alteration due to fitting, twisting or wearing of the insert can cause anomalous operation: in such case inserts must be systematically replaced.
After every procedure, unscrew the insert by means of the dynamometric wrench and sterilise.

USE OF INSERTS (FIG. 3)

- INSERT # AX-S Use: removal of supragingival tartar deposits on all tooth surfaces. For hard deposits, it is possible to set high power.
Power setting (min.→max): 3 ÷ 10
- INSERT # BX-S Use: removal of supragingival tartar, coatings and stains on all tooth surfaces.
Power setting (min.→max): 3 ÷ 10
- INSERT # PX-S Use: removal of subgingival tartar deposits on all tooth surfaces, including inter-proximal and sulcus areas.
Power setting (min.→max): 1 ÷ 6

The listed power values refer to settings which can be read on the potentiometer dial or LCD display of the power supply unit.
Suggested power settings are mean values given as a guide only.

REPLACEMENT OF O-RING

Should water leakage occur, replacement of the O-Ring (REF 001.18.03) on the handpiece connector is recommended.
With the aid of a pin (Fig.4), pull out the old O-ring and insert the new one on the corresponding groove.

PRELIMINARY PREPARATION

Maintenance processes are described according to ISO 17664 standard.
Before starting a process:

- Observe suitable personal protective measures against risk of infection.
- Remove all residual blood, cement or composite.
- Remove the insert from the instrument.
- Clean the medical device immediately after treatment of each patient.

Do not immerse the medical device into disinfectants solutions and do not use ultrasonic cleaning units.

CLEANING

- Manual external: Carefully clean the surfaces using a clean cloth soaked in a suitable substance. Use also a medium toothbrush and brush off under flowing hot tap water.
- Manual internal: Using a suitable spray product, follow corresponding instructions and spray the substance to the internal of the medical device. Immediately after internal cleaning, carry out the disinfection process or the drying process.
Use a thermosensitizer complying with ISO 15883-1 standard which is operated with alkaline cleaning agents at a pH of maximum 10.

DISINFECTION

- Manual external: Carefully clean the surfaces using a clean cloth soave in a suitable substance (i.e. O-phenylphenol or alcohols based on ethanol). Products not recommended: disinfectants containing benzalkonium chloride, acetone or glutaraldehyde.
- Manual internal: Using a suitable spray product, follow corresponding instructions and spray the substance to the internal of the medical device. Immediately after internal disinfection, carry out the sterilization process.
Use a thermosensitizer complying with ISO 15883-1 standard which is operated with alkaline disinfectant agents at a pH of maximum 10.

DRYING

- Manual: In general use dry and clean air complying to ISO 7494-2 standard.
- Automated: The drying phase is normally part of the cleaning program of the thermosensitizer.

STERILIZATION

Handpiece, wrench and inserts can be sterilised with fractioned initial vacuum phase, class B cycle according to the EN 13060 standard.
Sterilize at the nominal temperature of 134 °C / 2 bar (273 °F / 29 psi) for at least 3 min.
Use only demineralised or distilled water.
The handpiece must always be sterilised without the insert on. Use a sterilization bag.

The insert instead can be fit in and and sterilized together with the dynamometric wrench: this operation avoids any injury or contamination as it is not necessary to touch the insert.

The device is supplied "not sterile".
The device can be sterilised for at least 1000 times.

The processed device must be kept protected from dust and microbial contamination.

TRANSPORT AND STORAGE CONDITIONS

- Temperature: -10 ÷ 60 °C (14 ÷ 104 °F)
Relative humidity: 10 ÷ 85 %
Atmospheric pressure: 50 ÷ 106 KPa

SERVICE

In case of malfunction of the device and for any overhaul and repair work, please contact your usual supplier or apply directly to TKD.

GUARANTEE

The device carry a 24-month guarantee against all defects of construction.
Guarantee validity is subject to specific conditions.

REFUSE DISPOSAL

This device needs to be recycled. Electrical and electronic equipment may contain dangerous substances which constitute health and environmental hazards. The user must return the equipment to its dealer or establish direct contact with an approved body able to process and derive value from this type of equipment (European Directive 2002/96/EC).

TECHNICAL SPECIFICATIONS

- Standard conformity: ISO 18397, IEC 80601-2-60, CEI EN 60601-1, CEI EN 60601-1-2
- Operating frequency: 27000 ÷ 30000 Hz
- Classification: class 2 according to CEI EN 60601-1
- Protection against electrical shock: insert is a type B applied part according to IEC 80601-2-60
- Max temperature of applied part: 51 °C (124 °F)
- Operating Mode: Continuous (with spray)
Intermittent (without spray). Time ON: 2 min., Time OFF: 10 min.
- Weight without insert: 45 g.
Working temperature: 0 ÷ 40 °C (32 ÷ 104 °F)
Working humidity: 30 ÷ 90 %
External protection: Fuse T1AL
Input voltage: 24 Vac or 32 Vdc (USC60 electronic board)
Max. input current: 1 A
Max. output power: 9 W
Light source: LED 3500 °K, λ: 400 ÷ 600 nm

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